


Printing and/or Exporting a Falsified Medicines Report

The Falsified Medicines Report enables you to view the number of packs which returned an unknown product code, serial number and batch designation within your specified time period.

1. To begin, open the FMD History from the ProScript Connect Main Screen by clicking the **[ProScript Connect]** menu button and selecting the *FMD History* menu item.
2. From the FMD History, use the *Filters* to refine the grid to display the records to include in your report output.




The screenshot shows the ProScript Connect interface for FMD History. The top navigation bar includes the ProScript Connect logo, a search bar, and the user's name (Responsible Pharmacist Toby McElroy) and role (User Super Visor). The left sidebar contains various menu items like Redeem Owings, Reprints, Ordering, Overdue Scripts, MDS Prep, MDS Scripts Due, Stock Manager, MDS Manager, SDM, ETP, Script Q, and RMS. The main area displays a table of records with columns for Date Added, Date Sent, Pack De, GTIN, Serial Number, Batch Number, Action, and Status. The filters section shows a period of 'Last Year' from 01/01/2018 to 31/12/2018, with a status filter set to 'Show All'. The table contains 20 rows of data, with some rows showing 'DECOMMISSION' and 'SUCCESS' actions. A 'Result Description' box at the bottom indicates 'Unknown serial number.'


Date Added	Date Sent	Pack De	GTIN	Serial Number	Batch Number	Action	Status
28/11/2018 14:23:20	28/11/2018 14:23:21	[No Pac	00745119783481	PK158EDF33C3A432F47	TEST46AFC5E654	VERIFY	ERROR
28/11/2018 14:23:20	28/11/2018 14:23:20	[No Pac	00745119783481	PK0018EB2AA251D4749	TEST46AFC5E654	VERIFY	SUCCESS
28/11/2018 14:19:57	28/11/2018 14:19:59	[No Pac	00745119783481	PK1626698F7989E4340	TEST46AFC5E654	VERIFY	ERROR
13/11/2018 10:41:13	13/11/2018 10:41:15	Metfor	05012617009739	PK01099F0048485123T	TEST8133BCB979	VERIFY	ERROR
12/11/2018 15:36:31	12/11/2018 15:36:31	[No Pac	53106583345148	PK3261A6E454037A04D	TESTE7002D7EDF	VERIFY	ERROR
12/11/2018 15:34:17	12/11/2018 15:34:17	[No Pac	44714821587889	PK350D936A1A883FC45	TESTOD31C526D0	VERIFY	ERROR
12/11/2018 15:33:29	12/11/2018 15:33:29	[No Pac	44714821587889	PK3261A6E454037A04D	TESTOD31C526D0	VERIFY	ERROR
12/11/2018 15:32:31	12/11/2018 15:32:31	[No Pac	44714821587889	PK3261A6E454037A04D	TESTOD31C526D0	VERIFY	ERROR
12/11/2018 15:27:54	12/11/2018 15:27:54	[No Pac	44714821587889	PK0232A9CE5F57F8F47	TESTOD31C526D0	DECOMMISSION	SUCCESS
12/11/2018 15:19:04	12/11/2018 15:19:04	[No Pac	44714821587889	PK0200C449D8531444C	TESTOD31C526D0	DECOMMISSION	SUCCESS
12/11/2018 15:18:51	12/11/2018 15:18:51	[No Pac	44714821587889	PK022EBE93C8F9F8A46	TESTOD31C526D0	DECOMMISSION	ERROR
12/11/2018 14:45:06	12/11/2018 14:45:06	[No Pac	44714821587889	PK022EBE93C8F9F8A46	TESTOD31C526D0	VERIFY	ERROR
12/11/2018 14:44:57	12/11/2018 14:44:57	[No Pac	44714821587889	PK0200C449D8531444C	TESTOD31C526D0	VERIFY	SUCCESS
12/11/2018 14:25:49	12/11/2018 14:25:49	Atenolc	05017007019206	PK01099F0048485962O	TEST8133BCB979	VERIFY	ERROR
12/11/2018 14:11:53	12/11/2018 14:11:53	[No Pac	44714821587889	PK022EBE93C8F9F8A46	TESTOD31C526D0	VERIFY	ERROR
12/11/2018 14:11:53	12/11/2018 14:11:53	[No Pac	44714821587889	PK0200C449D8531444C	TESTOD31C526D0	VERIFY	SUCCESS
12/11/2018 13:55:34	12/11/2018 13:55:35	[No Pac	44714821587889	PK0200C449D8531444C	TESTOD31C526D0	VERIFY	SUCCESS
05/11/2018 10:48:50	05/11/2018 10:48:51	[No Pac	44714821587889	PK021772AF1D8883F4B	TESTOD31C526D0	DECOMMISSION	SUCCESS
05/11/2018 10:48:50	05/11/2018 10:48:50	[No Pac	44714821587889	PK019B6F568B4062945	TESTOD31C526D0	DECOMMISSION	SUCCESS

Result Description
Unknown serial number.

3. Select the [F11 – Reports] button, and select the *Print Falsified Medicines Report* menu item from the drop-down.



4. The Falsified Medicines Report preview opens. From here, you can:
 - Use the print icon to print the report
 - Use the floppy disk icon to save the report to your computer.



Revision History

Version Number	Date	Revision Details	Author(s)
1.0	7 th February 2019		Joanne Hibbert-Gorst

Contact us

For more information contact:
Telephone: 0344 209 2601