


# Printing and/or Exporting an FMD Main Report

The FMD Main Report enables you to view the number of pack scans completed – verification, decommission and recommission – within your specified time period.


1. To begin, open the FMD History from the ProScript Connect Main Screen by clicking the **[ProScript Connect]** menu button and selecting the *FMD History* menu item.
2. From the FMD History, use the *Filters* to refine the grid to display the records to include in your report output.




Date Added	Date Sent	Pack De:	GTIN	Serial Number	Batch Number	Action	Status
28/11/2018 14:23:20	28/11/2018 14:23:21	[No Pac	00745119783481	PK158EDF33C3A432F47	TEST46AFC5E654	VERIFY	ERROR
28/11/2018 14:23:20	28/11/2018 14:23:20	[No Pac	00745119783481	PK0018EB2AA251D4749	TEST46AFC5E654	VERIFY	SUCCESS
28/11/2018 14:19:57	28/11/2018 14:19:59	[No Pac	00745119783481	PK1626698F7989E4340	TEST46AFC5E654	VERIFY	ERROR
13/11/2018 10:41:13	13/11/2018 10:41:15	Metfor	05012617009739	PK01099F0048485123T	TEST8133BCB979	VERIFY	ERROR
12/11/2018 15:36:31	12/11/2018 15:36:31	[No Pac	53106583345148	PK3261A6E454037A04D	TESTE7002D7EDF	VERIFY	ERROR
12/11/2018 15:34:17	12/11/2018 15:34:17	[No Pac	44714821587889	PK350D936A1A883FC45	TEST0D31C526D0	VERIFY	ERROR
12/11/2018 15:33:29	12/11/2018 15:33:29	[No Pac	44714821587889	PK3261A6E454037A04D	TEST0D31C526D0	VERIFY	ERROR
12/11/2018 15:32:31	12/11/2018 15:32:31	[No Pac	44714821587889	PK3261A6E454037A04D	TEST0D31C526D0	VERIFY	ERROR
12/11/2018 15:27:54	12/11/2018 15:27:54	[No Pac	44714821587889	PK0232A9CE5F57F8F47	TEST0D31C526D0	DECOMMISSION	SUCCESS
12/11/2018 15:19:04	12/11/2018 15:19:04	[No Pac	44714821587889	PK0200C449D8531444C	TEST0D31C526D0	DECOMMISSION	SUCCESS
12/11/2018 15:18:51	12/11/2018 15:18:51	[No Pac	44714821587889	PK022EBE93C8F9F8A46	TEST0D31C526D0	DECOMMISSION	ERROR
12/11/2018 14:45:06	12/11/2018 14:45:06	[No Pac	44714821587889	PK022EBE93C8F9F8A46	TEST0D31C526D0	VERIFY	ERROR
12/11/2018 14:44:57	12/11/2018 14:44:57	[No Pac	44714821587889	PK0200C449D8531444C	TEST0D31C526D0	VERIFY	SUCCESS
12/11/2018 14:25:49	12/11/2018 14:25:49	Atenolc	05017007019206	PK01099F0048485962O	TEST8133BCB979	VERIFY	ERROR
12/11/2018 14:11:53	12/11/2018 14:11:53	[No Pac	44714821587889	PK022EBE93C8F9F8A46	TEST0D31C526D0	VERIFY	ERROR
12/11/2018 14:11:53	12/11/2018 14:11:53	[No Pac	44714821587889	PK0200C449D8531444C	TEST0D31C526D0	VERIFY	SUCCESS
12/11/2018 13:55:34	12/11/2018 13:55:35	[No Pac	44714821587889	PK0200C449D8531444C	TEST0D31C526D0	VERIFY	SUCCESS
05/11/2018 10:48:50	05/11/2018 10:48:51	[No Pac	44714821587889	PK021772AF1D8883F48	TEST0D31C526D0	DECOMMISSION	SUCCESS
05/11/2018 10:48:50	05/11/2018 10:48:50	[No Pac	44714821587889	PK01986F56884062945	TEST0D31C526D0	DECOMMISSION	SUCCESS

**Result Description**  
Unknown serial number.

- Select the **[F11 – Reports]** button, and select the *Print Main Report* menu item from the drop-down.



- The FMD Main Report preview opens. From here, you can:
  - Use the print icon to print the report
  - Use the floppy disk icon to save the report to your computer.



**FMD Main Report**

ALLIANCE PHARMACY  
 16 Beecham Road, Watford, Hertfordshire, WD18 8LP  
 From: 01/01/2018 To: 31/12/2018 Total Packs: 57

Date Added	Date Sent	Pack Description	GTIN	Serial No.	Batch No.	Expiry Date	Action	Status	User
28/11/2018 14:23	28/11/2018 14:23		00745119783481	PK158EDF33C3A 432F47	TEST46AFC5E65 4	28/11/2023	Verify	Error	Super Visor
28/11/2018 14:23	28/11/2018 14:23		00745119783481	PK001BEB2AA25 1D4749	TEST46AFC5E65 4	28/11/2023	Verify	Success	Super Visor
28/11/2018 14:19	28/11/2018 14:19		00745119783481	PK1626698F7989 E4340	TEST46AFC5E65 4	28/11/2023	Verify	Error	test test
13/11/2018 10:41	13/11/2018 10:41	Metformin 500mg tabs (28) ACTAVIS	05012617009739	PK01099F00484 B5123T	TEST8133BCB97 9	07/11/2023	Verify	Error	Super Visor
12/11/2018 15:36	12/11/2018 15:36		53106583345148	PK3261A6E4540 37A04D	TESTE7002D7E DF	03/10/2023	Verify	Error	Super Visor
12/11/2018 15:34	12/11/2018 15:34		44714821587889	PK350D936A1A8 83FC45	TEST0D31C526D 0	03/10/2023	Verify	Error	Super Visor
12/11/2018 15:33	12/11/2018 15:33		44714821587889	PK3261A6E4540 37A04D	TEST0D31C526D 0	03/10/2023	Verify	Error	Super Visor
12/11/2018 15:32	12/11/2018 15:32		44714821587889	PK3261A6E4540 37A04D	TEST0D31C526D 0	03/10/2023	Verify	Error	Super Visor
12/11/2018 15:27	12/11/2018 15:27		44714821587889	PK0232A9CE5F5 7F8F47	TEST0D31C526D 0	03/10/2023	Decommission	Success	Super Visor
12/11/2018 15:19	12/11/2018 15:19		44714821587889	PK0200C449D85 31444C	TEST0D31C526D 0	03/10/2023	Decommission	Success	Super Visor
12/11/2018 15:18	12/11/2018 15:18		44714821587889	PK022EBE93C8F 9F8A46	TEST0D31C526D 0	14/11/2018	Decommission	Error	Super Visor
12/11/2018 14:45	12/11/2018 14:45		44714821587889	PK022EBE93C8F 9F8A46	TEST0D31C526D 0	14/11/2018	Verify	Error	Super Visor
12/11/2018 14:44	12/11/2018 14:44		44714821587889	PK0200C449D85 31444C	TEST0D31C526D 0	03/10/2023	Verify	Success	Super Visor
12/11/2018 14:44	12/11/2018 14:44	Atenolol 25mg tabs (28) TEVA	05017007019206	PK01099F00484	TEST8133BCB97	07/11/2023	Verify	Error	Super Visor

## Revision History

Version Number	Date	Revision Details	Author(s)
1.0	7 <sup>th</sup> February 2019		Joanne Hibbert-Gorst

### Contact us

For more information contact:  
**Telephone: 0344 209 2601**

Not to be reproduced or copied without the consent of AAH Pharmaceuticals Limited 2019.