


Stock Control Configuration

There are several key application settings that can be amended to provide you with different methods of stock management within ProScript Connect.



Application Setting	Available Values	
Available Ordering Systems	Free-type a list of ordering systems to appear in the <i>Order System</i> drop-down menu in the Stock Record Card. Each order method must be separated by commas without any spacing. Available values include; Auto Auto Suppress Reminder Replace Don't Order	
After a stock count, users are able to select the [Stock Counted] Beliability Days & Estimate Days button within the item's Stock Record Card to update the Stock Status to Counted. The values set within the Reliability Days and		



	Estimate Days application settings determine the number of days needed to lapse before the Stock Status is updated.	
	 Reliability Days: Defines the number of days between the Counted and Estimate statuses. Estimate Days: Defines the number of days between the Estimate and Uncounted statuses. 	
Consolidate Stock Level	 YES; Consolidated stock mode will combine stock levels and usage figures for all drugs in the same family and category. NO; Non-consolidated stock mode will show stock levels and usage figures at item level. 	
Inhibit Stock Control	Allows you to disable all editable fields within the Stock Record Card to prevent manual amendments.	
Record Stock Adjustment Reason	Enable this setting to prompt users to provide an adjustment reason when manually amending stock levels. This information can then be seen on the Stock Adjustment Report tab in the Stock Manager.	

Revision History

Version Number	Date	Revision Details	Author(s)
1.0	19 August 2019		Joanne Hibbert-Gorst

Contact us

For more information contact: **Telephone: 0344 209 2601**

Not to be reproduced or copied without the consent of AAH Pharmaceuticals Limited 2019.