



Printing a Preparation or Pack Note

To begin, open the Drug Notes Manager from the ProScript Connect Main Screen by clicking the [ProScript Connect] menu button and selecting the *Drug Notes Manager* menu item.



2. From the Drug Notes Manager, select either the *Preparation* or *Pack Note* tab and use the available *Filters* to search for, locate and highlight the note you wish to print. Select the **[F3 – Print Note]** button.





3. The note label print preview displays. Click the printer icon to print.

Revision History

Version Number	Date	Revision Details	Author(s)
1.0	22 March 2023		Joanne Hibbert-Gorst

Not to be reproduced or copied without the consent of AAH Pharmaceuticals Limited 2021.

Contact us

For more information contact: **Telephone: 0344 209 2601**