


Deleting/Restoring a Setting Template

Deleting a Setting Template

1. To begin, open the Setting Templates Manager from the ProScript Connect Main Screen by clicking the **[ProScript Connect]** menu button and selecting the *Setting Templates Manager* menu item.
2. From the Setting Templates Manager, highlight a template and select the **[F2 – Edit Selected Template]** button. The Edit Template Details form displays.



Edit Template Details - England Branch Settings

Template Details


Template Name *

Settings + F6 - Template Settings

Category	Name	Value
Default	IdleLockSeconds	150
Default	ShowRxSummary	NO

Template Deleted | Last edited on 07/06/2016 at 11:22 ✓ F10 - Save ✗ Cancel

3. Select the *Template Deleted* tick box. A DELETED watermark displays across the window.



Edit Template Details - England Branch Settings

Template Details

Template Name *

Settings + F6 - Template Settings


Category	Name	Value
Default	IdleLockSeconds	150
Default	ShowRxSummary	NO

Template Deleted | Last edited on 07/06/2016 at 11:22 ✓ F10 - Save ✗ Cancel

4. Select the **[F10 – Save]** button.

Restoring a Setting Template

1. To begin, open the Setting Templates Manager from the ProScript Connect Main Screen by clicking the **[ProScript Connect]** menu button and selecting the *Setting Templates Manager* menu item.
2. From the Setting Templates Manager, select the *Show Deleted* tick box. Any deleted setting templates now display with a **DELETED** tag.



3. Highlight the template you wish to restore and select the **[F2 – Edit Selected Template]** button. The Edit Template Details window opens with a *DELETED* watermark.
4. Deselect the *Template Deleted* tick box in the bottom left-hand corner and select the **[F10 – Save]** button.

Edit Template Details - BNF Settings
— □ ×

Template Details

Template Name *

Settings + F6 - Template Settings

Category	Name	Value
BNF	BNFEnabled	NO

Template Deleted | Last edited on 01/07/2020 at 12:00
✓ F10 - Save
✗ Cancel

Revision History

Version Number	Date	Revision Details	Author(s)
1.0	15/08/2023		Joanne Hibbert-Gorst

Not to be reproduced or copied without the consent of AAH Pharmaceuticals Limited 2023.

Contact us

For more information contact:
Telephone: 0344 209 2601